

SUPPLEMENTAL MATERIAL


Salimi et al., <http://www.jem.org/cgi/content/full/jem.20130351/DC1>

Figure S1. Gating strategies for the identification of human ILC2 in blood and skin using flow cytometry. (A) Gating strategy for sorting of ILC2 from the blood. Total PBMCs were stained for CD45, lineage markers, IL-7R α , and CRTH2, and ILC2s were gated as CD45⁺Lin⁻IL-7R α ⁺CRTH2⁺. Data are representative of five independent experiments ($n = 18$ donors). (B) Cellular infiltrate into the skin after HDM challenge in allergic and nonallergic donors. Shown is a gating strategy to show different populations of cells that infiltrated human skin blisters after HDM allergen or PBS challenge.